

# Immunoblots of Markers of the Unfolded Protein Response: Method for Preparing Cytoplasmic and Nuclear Extracts

From Dr. Heather Harding of the laboratory of Dr. David Ron

## Solutions:

### Harvest buffer:

10 mM HEPES pH 7.9  
50 mM NaCl  
0.5 M Sucrose  
0.1 mM EDTA  
0.5% Triton-X100

### For 50ml:

0.5 ml 1 M HEPES  
0.5 ml 5 M NaCl  
8.56 g sucrose  
10  $\mu$ l 0.5 M EDTA  
2.5 ml 10% Triton-X100

\*freshly add to harvest buffer just before use:

1 mM DTT  
10 mM tetrasodium pyrophosphate  
100 mM NaF  
17.5 mM  $\beta$ -glycerophosphate  
1 mM PMSF  
4  $\mu$ g/ml Aprotinin  
2  $\mu$ g/ml Pepstatin A

### Buffer A:

10 mM HEPES pH 7.9  
10 mM KCL  
0.1 mM EDTA  
0.1 mM EGTA

### For 50ml:

0.5 ml 1 M HEPES  
250  $\mu$ l 2 M KCL  
10  $\mu$ l 0.5 M EDTA  
10  $\mu$ l 0.5 M EGTA

Freshly add to buffer A just before use:

1 mM DTT  
1 mM PMSF  
4  $\mu$ g/ml Aprotinin  
2  $\mu$ g/ml Pepstatin A

### Buffer C:

10 mM HEPES pH 7.9  
500 mM NaCl  
0.1 mM EDTA  
0.1 mM EGTA  
0.1% IGEPAL(NP40)

### For 50ml 1X:

0.5ml 1 M HEPES  
5 ml 5 M NaCl  
10  $\mu$ l 0.5 M EDTA  
10  $\mu$ l 0.5 M EGTA  
500  $\mu$ l 10% IGEPAL

### for 50ml 2X

1 ml 1 M HEPES  
10 ml 5 M NaCl  
20  $\mu$ l 0.5 M EDTA  
20  $\mu$ l 0.5 M EGTA  
1 ml 10% IGEPAL

Freshly add to buffer C just before use:

1 mM DTT  
1 mM PMSF  
4  $\mu$ g/ml Aprotinin  
2  $\mu$ g/ml Pepstatin A

Ice cold PBS and PBS-1mM EDTA

## Method:

- Place plates on ice.
- Wash 2X with cold PBS.
- Add 1 ml PBS-EDTA and scrape cells. Transfer to 1.5ml tube.
- Pellet at 3000 rpm for 5 min.
- Resuspend in 250-500  $\mu$ l (at least 6 vol) of Harvest buffer.
- Incubate on ice for 5 min.
- Pellet at 1000 rpm in a table top swinging bucket rotor (for example by resting microfuge tubes in 15ml tubes while spinning) for 10 min to pellet nuclei @ 4°C (preferable). Alternatively spin in a fixed angle rotor at 2000 rpm. Note that the nuclei sometimes smear up the side of the tube in the fixed angle rotor.
- Transfer supernate to new tube. For best results, clear at 14,000 rpm for 15 min. and transfer supernate to new tube = "**Cytoplasmic/membrane extract**".  
Save for PERK or IRE IP or eIF2 $\alpha$  western.
- Wash/resuspend nuclear pellet in 500  $\mu$ l buffer A (optional).
- Pellet at 1000 rpm in swinging bucket rotor 5 min 4°C.
- Remove and discard sup.
- Add 4X vol buffer C. (*For a more concentrated extract use 2 vol buffer C made with 1M NaCl.*)
- Vortex 15 min at 4°C. (in cold room) Start on high speed vortex to loosen pellet, then turn to medium.
- Pellet at 14,000 rpm 10 min 4°C.
- Transfer sup to new tube = "**Nuclear extract**".

If nuclear extract will be used for functional assays such as gel shift (EMSA), then 10% glycerol should be added prior to freezing.

Measure protein concentrations of cytoplasmic and nuclear extracts.

Load ~40-50  $\mu$ g of nuclear extract for XBP1, ATF4 and CHOP Western (20  $\mu$ g is enough for many cell types, but loading more is better.)

Load 25-40  $\mu$ g of cytoplasmic extract for Phospho-eIF2 $\alpha$  Western.

Use 500-1500  $\mu$ g of cytoplasmic extract for PERK or IRE immunoprecipitation followed by Western (the more the better).

## Antibodies for Immunoblots:

Antibody	MW of protein	Source	Dilution
CHOP	29 kD	rabbit	1:5000
XBP-1	54 kD	rabbit	1:2000
PERK	125 kD*	rabbit	1:1000-3000
IRE1	110 kD*	rabbit	1:1000-3000
ATF4	42 kD	rabbit	1:3000

\*phosphorylated forms migrate slightly slower on gel